

# Biosafety Products: Sharps and Histology Containers

Sharps Container



Medical Waste Can



Histology Container



Biohazard Bags



# Biosafety Products: Sharps and Histology Containers



## Sharps Container

- Made of heavy duty plastic - puncture resistant.
- Properly labelled to warn of hazardous waste inside.
- Has a fill line.
- Leak-proof on the sides and bottom.

# Biosafety Products: Sharps and Histology Containers



## Medical Waste Can

- A color coded disposal unit where some infectious waste can be temporarily dumped.
- Must be lined with a yellow or red medical waste bag.

# Biosafety Products: Sharps and Histology Containers



## Biohazard Bag

- The two pouch design separates the document and the specimen to avoid contamination in case of leak.
- Available in reclosable zip top and self adhesive tamper proof strip.

# Biosafety Products: Sharps and Histology Containers



## Histology Containers

- Polycarbonate material that can withstand the corrosive chemicals and is suitable for the transport of specimens.
- Airtight and leak proof to prevent spillage of formalin.