



Trusted
Reliable
Unparelled

TruLab Hub Cutters

Safety Box for Disposal of
Hypodermic
Needles



TruLaboratories Corporation is promoting Injection Safety in the Philippine market.

TruLab Hub Cutters

Simple • Convenient • User Friendly

The TruLab Hub Cutters are portable, disposable, point-of-use needle hub destruction devices. The devices are designed to cut the syringe needle at the hub and enclose the needle securely in the hard plastic container. The Hub Cutters are disposed of as a whole unit when full. **A simple, convenient and affordable solution for medical sharps waste disposal in emerging markets.**

Product Benefits and Features

Benefits	Features
Safety	Cuts needle at the hub - no needle sharps remaining on the syringe after the cut
	Separates needle and syringe at the point-of-use
	Minimal splatter
	Puncture-resistant container
	Encloses needles securely in the hard plastic container
	Designated fill line
	Yellow color - color coding for biohazardous waste
	Container cap can be closed between uses
	Container cap can be permanently locked prior to disposal of product
	Disposal of entire container avoids contact with sharps
Convenience	Does not require electricity to operate
	No assembly required before or after use
	No disassembly required prior to disposal
User friendly	Compatible with a broad range of syringes
	Large opening and surrounding wall facilitate needle insertion
	Multilingual labeling
Outreach use	Portable, point-of-use, compact and lightweight
	Easy to use and carry around
Reduces the volume of medical sharps waste	Cuts and holds 400-600 needles, depending on the needle size.
Easy to dispose	Can be incinerated or buried in pit
Accessible	Affordable

